

BACTERIA TESTS

DIPSLIDES, TUBE TESTS, INCUBATORS

EASY AND RELIABLE

FLEXIBLE PADDLE

VERSIONS

TTC (TOTAL VIABLE COUNT)

YEAST-MOULD-FUNGI/TTC

COLIFORMS/TTC

E. COLI - PSEUD. AE./TTC

PSEUDOMON. SPEC./TTC

SULPHATE RED. BACT.

NITRITE RED. BACT./TTC

DIPSLIDES FOR VARIOUS SORTS AND COMBINATIONS OF BACTERIA



BACTEDT057
INCUBATOR FOR 2 DIPSLIDES



BACTEDT055
INCUBATOR FOR 8 DIPSLIDES



BACTEDT059
INCUBATOR FOR 10 DIPSLIDES



BACTEDT052
INCUBATOR FOR 25 DIPSLIDES



BACTEDT093
INCUBATOR FOR 45 DIPSLIDES



WATER-I.D.®

WATER TESTING EQUIPMENT ●●●

BACTERIA TESTS

DIPSLIDES, TUBE TESTS, INCUBATORS

FEATURES

The definition of a dipslide is a means of testing the microbial content of liquids or surfaces. It consists of a plastic carrier bearing a sterile culture medium which can be dipped in the liquid or pressed on to the surface to be sampled. It is then incubated for 48 hours at a specific temperature to allow microbial growth.

The resulting microbial colonies are estimated by reference to a chart. If the treatment programme is effective, the dipslide count will be consistently low (<10³). If a high result is obtained, the test should be repeated and the treatment programme checked.

Water-i.d.® dipslides for bacterial testing are provided as 10 dipslides in a box.

Water-i.d.® offers a variety of incubators along with the dipslide range.

Water-i.d.® Dipslides for Bacteria testing



Example:
Incubator BactEDT055
for 8 dipslides

DIPSLIDES

Each Water-i.d.® Dipslide pack contains 10 dipslides along with a comparison chart to obtain the test result.

Item code	Product Description
BactD001	TTC/TTC Dipslide Tests for Total Viable Count in waters and surfaces. Flexible paddle for effective surface contact.
BactD002	Malt/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count in waters and surfaces. Flexible paddle for effective surface contact.
BactD003	RBS/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count for food and manufacturing industry. Flexible paddle for effective surface contact.
BactD004	MAG/TTC Dipslide Tests for Coliforms and Total Viable count in process water. Flexible paddle for effective surface contact.
BactD005	Chromogenic/TTC Dipslide Tests for E.coli, Pseudomonas Aeruginosa and Total Viable Count for enumeration of E.coli. Flexible paddle for effective surface contact.
BactD006	PDMMAC Dipslide Tests for Pseudomonas species and Coliforms for potable water, pools and spas. Flexible paddle for effective surface contact.
BactD007	PDW/TTC Dipslide Tests for Pseudomonas species and Total Viable Count for Closed water systems, pools and spas. Flexible paddle for effective surface contact.
BactD008	SRB Tube Tests for Sulphate reducing bacteria in Closed water systems. Flexible paddle for effective surface contact.
BactD009	NRB Tube Tests for Nitrite reducing bacteria in Closed water systems. Flexible paddle for effective surface contact.

INCUBATORS

Item code	Product Description	Voltage
BactEDT057	Incubator for 2 Dipslide Tests. Preset to 30°C.	110 / 220 V
BactEDT055	Incubator for 8 Dipslide Tests. With temperature gauge and adjustable temperature.	12 / 110 / 220 V
BactEDT059	Incubator for 10 Dipslide Tests. With temperature gauge and adjustable temperature. Ideal for NRB/SRB tests. Only offered with 240 V.	12 / 110 / 220 V
BactEDT052	Incubator for 25 Dipslide Tests. With temperature gauge and adjustable temperature.	12 / 110 / 220 V
BactEDT093	Incubator for 45 Dipslide Tests. With temperature gauge and adjustable temperature.	12 / 110 / 220 V

We gladly Private Label our products for you!
Write to sales@water-id.com to receive your individual quotation!

WATER-I.D.®

WATER TESTING EQUIPMENT ●●●

BACTERIA TESTS

DIPSLIDES, TUBE TESTS, INCUBATORS



DIPSLIDE PACKS

Dipslide Packs are always 10 double sided dipslides with flexible paddle for effective surface contact.

Package includes comparison chart and instructions.

Item No. Product Description

BactD001	TTC/TTC Dipslide Tests for Total Viable Count in waters and surfaces.
BactD002	Malt/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count in waters and surfaces.
BactD003	RBS/TTC Dipslide Tests for Yeasts, Moulds, Fungi and Total Viable Count for food and manufacturing industry.
BactD004	MAC/TTC Dipslide Tests for Coliforms and Total Viable count in process water.
BactD005	Chromogenic/TTC Dipslide Tests for E.coli, Pseudomonas Aeruginosa and Total Viable Count for enumeration of E.coli.
BactD006	PDM/MAC Dipslide Tests for Pseudomonas species and Coliforms for potable water, pools and spas.
BactD007	PDM/TTC Dipslide Tests for Pseudomonas species and Total Viable Count for Closed water systems, pools and spas.
BactD008	SRB Tube Tests for Sulphate reducing bacteria in Closed water systems.
BactD009	NRB Tube Tests for Nitrite reducing bacteria in Closed water systems.



INCUBATORS

Item No. Product Description

Item No.	Product Description	Maximum dipslide qty	Voltage
BactEDT057	Incubator for 2 Dipslide Tests. Preset to 30°C.	2	110 / 220 V
BactEDT055	Incubator for 8 Dipslide Tests. With temperature gauge and adjustable temperature.	8	12 / 110 / 220 V
BactEDT059	Incubator for 10 Dipslide Tests. With temperature gauge and adjustable temperature. Ideal for NRB/SRB tests. Only offered with 240 V.	10	12 / 110 / 220 V
BactEDT052	Incubator for 25 Dipslide Tests. With temperature gauge and adjustable temperature.	25	12 / 110 / 220 V
BactEDT093	Incubator for 45 Dipslide Tests. With temperature gauge and adjustable temperature.	45	12 / 110 / 220 V
BactEDT-12V	12V car-lead for BactEDT057/BactEDT055/ BactEDT052/BactEDT093		12 V

We gladly Private Label our products for you!
Write to sales@water-id.com to receive your individual quotation!